



Anti-PTPN23 polyclonal antibody (CABT-BL3059)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide of Human HDPTP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic
Format	Liquid
Size	100 µg
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Tyrosine-protein phosphatase non-receptor type 23 is an enzyme that in humans is encoded by

the PTPN23 gene.

GENE INFORMATION

Entrez Gene ID	25930
----------------	-----------------------

Protein Refseq	NP_056281
----------------	---------------------------

UniProt ID	Q9H3S7
------------	------------------------
