



Anti-PUS3 polyclonal antibody (CABT-BL3065)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human PUS3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The protein encoded by this gene catalyzes the formation of tRNA pseudouridine from tRNA uridine at position 39 in the anticodon stem and loop of transfer RNAs. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2012]
--------------	--

GENE INFORMATION

Entrez Gene ID [83480](#)

Protein Refseq [NP_112597.3](#)

UniProt ID [A0A087WY59](#)
