



## Anti-PUT (aa 128-147) polyclonal antibody (CABT-BL3066)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                    |                                                                                                                                  |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Immunogen          | Synthetic peptide: KSTTEATTQVPKEKTQDGSN, corresponding to C terminal amino acids 128-147 of <i>Drosophila melanogaster</i> Punt. |
| Isotype            | IgG                                                                                                                              |
| Source/Host        | Rabbit                                                                                                                           |
| Species Reactivity | <i>Drosophila</i>                                                                                                                |
| Purification       | Immunogen affinity purified                                                                                                      |
| Conjugate          | Unconjugated                                                                                                                     |
| Applications       | IP, WB, ICC                                                                                                                      |
| Format             | Liquid                                                                                                                           |
| Preservative       | 0.02% Sodium Azide                                                                                                               |
| Storage            | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.                                            |

### GENE INFORMATION

|                |                           |
|----------------|---------------------------|
| Entrez Gene ID | <a href="#">41772</a>     |
| Protein Refseq | <a href="#">NP_731926</a> |
| UniProt ID     | <a href="#">Q24468</a>    |