



Anti-PXMP2 (aa 47-76) polyclonal antibody (CABT-BL3067)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 47-76 of Human PXMP2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Peroxisome membrane.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	5827
Protein Refseq	NP_061133

UniProt ID

[Q9NR77](#)
