



Anti-PXMP4 polyclonal antibody (CABT-BL3068)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthetic peptide selected from the central region of Human PXMP4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Peroxisome membrane; Multi pass membrane protein.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	Peroxisomal membrane protein 4 is a protein that in humans is encoded by the PXMP4 gene.
--------------	--

GENE INFORMATION

Entrez Gene ID [11264](#)

Protein Refseq [NP_899634](#)

UniProt ID [Q9Y6I8](#)
