



Anti-PXMP4 (N-terminal) polyclonal antibody (CABT-BL3069)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthesized peptide derived from the N-terminal region of human PXMP4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Peroxisome membrane; Multi pass membrane protein.
Format	Liquid
Buffer	50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Peroxisomal membrane protein 4 is a protein that in humans is encoded by the PXMP4 gene.
--------------	--

GENE INFORMATION

Entrez Gene ID [11264](#)

Protein Refseq [NP_899634](#)

UniProt ID [Q9Y6I8](#)
