



Anti-QSOX1 (internal region) polyclonal antibody (CABT-BL3074)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from an internal domain of Mouse Quiescin Q6
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-P, WB
Cellular Localization	Golgi apparatus membrane and Secreted; extracellular space. Found in the extracellular medium of quiescent cells but is not found in proliferating cells.
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Sulfhydryl oxidase 1 is an enzyme that in humans is encoded by the QSOX1 gene.
---------------------	--

GENE INFORMATION

Entrez Gene ID [104009](#)

Protein Refseq [NP_001020116](#)

UniProt ID [Q8BND5](#)
