



# Anti-QSOX1 (aa 231-587) polyclonal antibody (CABT-BL3075)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant fragment, corresponding to a region within amino acids 231-587 of Human Quiescin Q6.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Golgi apparatus membrane and Secreted; extracellular space. Found in the extracellular medium of quiescent cells but is not found in proliferating cells.
<b>Format</b>	Liquid
<b>Buffer</b>	10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction** Sulphydryl oxidase 1 is an enzyme that in humans is encoded by the QSOX1 gene.

# GENE INFORMATION

Entrez Gene ID [5768](#)

---

Protein Refseq [NP\\_001004128](#)

---

UniProt ID [O00391](#)

---