



Anti-R3HDM1 polyclonal antibody (CABT-BL3078)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human R3HDM1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	R3HDM1 (R3H domain containing 1) is a protein-coding gene. GO annotations related to this gene include nucleic acid binding. An important paralog of this gene is R3HDM2.
--------------	---

GENE INFORMATION

Entrez Gene ID [23518](#)

Protein Refseq [NP_056176.2](#)

UniProt ID [Q15032](#)
