



Anti-RAB31 (aa 128-177) polyclonal antibody (CABT-BL3080)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C-terminal amino acids 128-177 (EVPLKDAKEY AESIGAIVVE TSAKNAINIE ELFQGISRQI PPLDPHENGN) of Human RAB31.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane; Lipid-anchor; Cytoplasmic side.
Format	Liquid
Buffer	98% PBS, 2% Sucrose
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Small GTP-binding proteins of the RAB family, such as RAB31, play essential roles in vesicle and granule targeting (Bao;et al., 2002 (PubMed 11784320)).
---------------------	--

GENE INFORMATION

Entrez Gene ID [11031](#)

Protein Refseq [NP_006859](#)

UniProt ID [Q13636](#)
