



Anti-RAB31 (full length) polyclonal antibody (CABT-BL3081)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length Human RAB31 protein (AAH01148.1).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane; Lipid-anchor; Cytoplasmic side.
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Small GTP-binding proteins of the RAB family, such as RAB31, play essential roles in vesicle and granule targeting (Bao;et al., 2002 (PubMed 11784320)).
--------------	--

GENE INFORMATION

Entrez Gene ID [11031](#)

Protein Refseq [NP_006859](#)

UniProt ID [Q13636](#)
