



# Anti-RAB31 (full length) polyclonal antibody (CABT-BL3082)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Full length protein, corresponding to amino acids 1-195 of Human RAB31 (AAH01148.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane; Lipid-anchor; Cytoplasmic side.
Format	Liquid
Buffer	1X PBS, pH 7.2
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Small GTP-binding proteins of the RAB family, such as RAB31, play essential roles in vesicle and granule targeting (Bao;et al., 2002 (PubMed 11784320)).
--------------	--

# GENE INFORMATION

Entrez Gene ID [11031](#)

---

Protein Refseq [NP\\_006859](#)

---

UniProt ID [Q13636](#)

---