



# Anti-RAB39 (full length) polyclonal antibody (CABT-BL3084)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                       |                                                                                       |
|-----------------------|---------------------------------------------------------------------------------------|
| Immunogen             | Full length Human RAB39 protein                                                       |
| Isotype               | IgG                                                                                   |
| Source/Host           | Mouse                                                                                 |
| Species Reactivity    | Human                                                                                 |
| Purification          | Whole antiserum                                                                       |
| Conjugate             | Unconjugated                                                                          |
| Applications          | WB                                                                                    |
| Cellular Localization | Cell membrane.                                                                        |
| Format                | Liquid                                                                                |
| Buffer                | Whole serum                                                                           |
| Preservative          | None                                                                                  |
| Storage               | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

## BACKGROUND

|              |                                    |
|--------------|------------------------------------|
| Introduction | RAB39A, member RAS oncogene family |
|--------------|------------------------------------|

## GENE INFORMATION

Entrez Gene ID [54734](#)

---

Protein Refseq [NP\\_059986](#)

---

UniProt ID [Q14964](#)

---