



# Anti-RAB39 (aa 1-100) polyclonal antibody (CABT-BL3085)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH derived from within residues 1 - 100 of Mouse RAB39.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane.
Format	Liquid
Buffer	PBS pH: 7.40
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

Introduction	RAB39A, member RAS oncogene family
--------------	------------------------------------

# GENE INFORMATION

Entrez Gene ID [270160](#)

---

Protein Refseq [NP\\_780771](#)

---

UniProt ID [Q0PD15](#)

---