



Anti-RAB6C (internal region) polyclonal antibody (CABT-BL3087)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from an internal region of human RAB6C.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

BACKGROUND

Introduction	Ras-related protein Rab-6C is a protein that in humans is encoded by the RAB6C gene.
--------------	--

GENE INFORMATION

Entrez Gene ID [84084](#)

Protein Refseq [NP_115520](#)

UniProt ID [Q9H0N0](#)
