



Anti-RAB6C (aa 118-167) polyclonal antibody (CABT-BL3088)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------------|---|
| Immunogen | Synthetic peptide corresponding to a region within internal amino acids 118-167 (DVIITLVGNRTDLADKRQVSVEEGERKAKGLNVTFIETR AKAGYNVKQLF) of Human RAB6C. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Immunogen affinity purified |
| Conjugate | Unconjugated |
| Applications | WB |
| Format | Liquid |
| Concentration | Lot specific |
| Buffer | 98% PBS, 2% Sucrose |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |

BACKGROUND

| | |
|--------------|--|
| Introduction | Ras-related protein Rab-6C is a protein that in humans is encoded by the RAB6C gene. |
|--------------|--|

GENE INFORMATION

Entrez Gene ID [84084](#)

Protein Refseq [NP_115520](#)

UniProt ID [Q9H0N0](#)
