



Anti-RABA4B (aa 201-219) polyclonal antibody (CABT-BL3089)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: GKKILIPGSGQEIPAKTST conjugated to KLH, corresponding to C terminal amino acids 201-219 of Arabidopsis thaliana AtRabA4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana, Soybean
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	0.01M PBS, pH 7.4
Preservative	15mM Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	GTP-binding protein ATGB3.
--------------	----------------------------

GENE INFORMATION

Entrez Gene ID [830160](#)

Protein Refseq [NP_195709](#)

UniProt ID [Q9SMQ6](#)
