



Anti-RAD54B (aa 86-186) polyclonal antibody (CABT-BL3094)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant tagged fusion protein: YCFEAPTLAT LDPPTHVHSA PKEVAVSKEQ EEKSDSLVKY FSVVWCKPSK KKHKKWEGDA VLIVKGKSF I LKNLEGKDIG RGIGYKFKELEKIEEGQTLM, corresponding to amino acids 86-186 of Human RAD54B
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	The RAD52 epistasis group genes are involved in homologous recombination, and they are conserved from yeast to humans. RAD54B is involved in DNA repair and mitotic recombination.
---------------------	--

GENE INFORMATION

Entrez Gene ID [25788](#)

Protein Refseq [NP_001192191](#)

UniProt ID [O95073](#)
