



Anti-RAD54B (aa 861-910) polyclonal antibody (CABT-BL3095)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the internal amino acids 861-910 (NSLKPLSMSQ LKQWKHFSGD HLNLTDPFLE RITENVSFIF QNITTQATGT) of human RAD54B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, ICC/IF
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction The RAD52 epistasis group genes are involved in homologous recombination, and they are

conserved from yeast to humans. RAD54B is involved in DNA repair and mitotic recombination.

GENE INFORMATION

Entrez Gene ID	25788
----------------	-----------------------

Protein Refseq	NP_001192191
----------------	------------------------------

UniProt ID	O95073
------------	------------------------
