



Anti-RAPGEF6 (aa 625-675) polyclonal antibody (CABT-BL3098)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between residues 625 and 675 of human RAGEF2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Cellular Localization	Cell Membrane and Cytoplasmic. Note: Upon binding to M-Ras, it translocates to the plasma membrane.
Format	Liquid
Buffer	0.1% BSA, Tris buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	RAPGEF6 (Rap guanine nucleotide exchange factor (GEF) 6) is a protein-coding gene. Diseases associated with RAPGEF6 include schizophrenia. GO annotations related to this gene include Ras GTPase binding and guanyl-nucleotide exchange factor activity. An important
---------------------	--

paralog of this gene is RAPGEFL1.

GENE INFORMATION

Entrez Gene ID	51735
----------------	-----------------------

Protein Refseq	NP_001157860
----------------	------------------------------

UniProt ID	B2RTU6
------------	------------------------
