



Anti-RBM18 (aa 101-150) polyclonal antibody (CABT-BL3106)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 101-150 (VRWAHAQVKR YDHNKNDKIL PISLEPSSST EPTQSNLSVT AKIKAIEAKL) of Human RBM18.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Concentration	Lot specific
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [92400](#)

Protein Refseq [NP_149108](#)

UniProt ID [Q96H35](#)
