



Anti-RBM26 polyclonal antibody (CABT-BL3107)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human RBM26.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	RBM26 (RNA Binding Motif Protein 26) is a Protein Coding gene. GO annotations related to this gene include RNA binding and nucleotide binding. An important paralog of this gene is RBM27.
---------------------	--

GENE INFORMATION

Entrez Gene ID [64062](#)

Protein Refseq [NP_071401](#)

UniProt ID [A0A087X0H9](#)
