



# Anti-RCOR3 (aa 1-50) polyclonal antibody (CABT-BL3118)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Immunogen** A region within synthetic peptide: MRVGAEQAR IPEFDPGATK YTDKDNGGML VWSPYHSIPD AKLDEYIAIA, corresponding to N terminal amino acids 1-50 of Human CoREST3 / RCOR3

**Isotype** IgG

**Source/Host** Rabbit

**Species Reactivity** Human

**Purification** Protein A purified

**Conjugate** Unconjugated

**Applications** WB, IHC-P, ELISA

**Cellular Localization** Nuclear

**Format** Liquid

**Buffer** 2% Sucrose, PBS

**Preservative** None

**Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

**Entrez Gene ID** [55758](#)

Protein Refseq [NP\\_001129697](#)

---

UniProt ID [Q9P2K3](#)

---