



Anti-Reg1B polyclonal antibody [Biotin] (CABT-BL3126)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Conjugate	Biotin
Applications	WB
Preservative	None
Storage	-20°C, Avoid freeze/thaw

BACKGROUND

Introduction	<p>This gene is a type I subclass member of the Reg gene family. The Reg gene family is a multigene family grouped into four subclasses, types I, II, III and IV based on the primary structures of the encoded proteins. This gene encodes a protein secreted by the exocrine pancreas that is highly similar to the REG1A protein. The related REG1A protein is associated with islet cell regeneration and diabetogenesis, and may be involved in pancreatic lithogenesis. Reg family members REG1A, REGL, PAP and this gene are tandemly clustered on chromosome 2p12 and may have arisen from the same ancestral gene by gene duplication.</p>
--------------	---

GENE INFORMATION

Entrez Gene ID	5968
----------------	----------------------

Protein Refseq [NP_006498](#)

UniProt ID [P48304](#)
