



Anti-REG4 polyclonal antibody (CABT-BL3129)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC, WB
Preservative	None
Storage	-20°C, Avoid freeze/thaw

BACKGROUND

Introduction	REG protein was shown to be stimulated during the regeneration of pancreatic islets. Since then, many Reg-related proteins have been identified in humans and other animals. In human, the four REG family genes, i.e., REG 1 alpha, REG 1 beta, REG-related s
--------------	--

GENE INFORMATION

Entrez Gene ID	83998
Protein Refseq	NP_001152824
UniProt ID	A0A024R0M1