



Anti-REI1 (C-terminal) polyclonal antibody (CABT-BL3130)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from the C terminal domain of Pre 60S factor REI1 from <i>Saccharomyces cerevisiae</i> . (Bakers yeast).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ICC, WB
Cellular Localization	Cytoplasmic
Format	Liquid
Buffer	Preservative: None
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [852569](#)

Protein Refseq [NP_009825](#)

UniProt ID [P38344](#)
