



Anti-REM2 (full length) polyclonal antibody (CABT-BL3132)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length Human REM2 protein (AAH35663, 1 a.a. ~ 331 a.a).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, ICC/IF
Cellular Localization	Cell Membrane
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	161253
Protein Refseq	NP_775798

UniProt ID

[Q8LYK8](#)
