



Anti-REP15 polyclonal antibody (CABT-BL3133)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human REP15.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	REP15 is a binding partner of the RAB GTPase family member RAB15 that facilitates transferrin receptor (TFRC; MIM;190010) recycling from the endocytic recycling compartment (Strick and Elferink, 2005 (PubMed 16195351)).
--------------	---

GENE INFORMATION

Entrez Gene ID [387849](#)

Protein Refseq [NP_001025045](#)

UniProt ID [Q6BDI9](#)
