



Anti-RETNLA (full length) polyclonal antibody (CABT-BL3135)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant full length protein (Mouse)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr, ELISA, WB
Cellular Localization	Secreted.
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Probable hormone.
---------------------	-------------------

GENE INFORMATION

Entrez Gene ID [57262](#)

Protein Refseq [NP_065255](#)

UniProt ID [G3X8W2](#)
