



Anti-RETNLA (aa 31-44) polyclonal antibody (CABT-BL3136)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to amino acids 31-44 (VENKVKELLANPAN) of Mouse RELM alpha.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted.
Format	Liquid
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Probable hormone.
---------------------	-------------------

GENE INFORMATION

Entrez Gene ID [57262](#)

Protein Refseq [NP_065255](#)

UniProt ID [G3X8W2](#)
