



Anti-RETNLA (full length) polyclonal antibody [Biotin] (CABT-BL3137)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant full length protein (Mouse)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Biotin
Applications	FC, ELISA, WB
Reconstitution	Reconstitute to a concentration of 50µg/ml with sterile PBS solution containing 0.1% BSA.
Cellular Localization	Secreted.
Format	Lyophilised
Size	50 µg
Buffer	0.1% BSA, PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

GENE INFORMATION

Entrez Gene ID	57262
----------------	-----------------------

Protein Refseq	NP_065255
----------------	---------------------------

UniProt ID	G3X8W2
------------	------------------------
