



Anti-REV3L (aa 418-518) polyclonal antibody (CABT-BL3138)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fusion protein: LAGLESDGYR GERNRMPSPC RSFGNNKYPQ NSDDEENEPQ IEKEEMELSL VMSQRWDSNI EEHCAKKRSL CRNTHRSSTE DDDSSSGEEM EWSDNSLLLA, corresponding to amino acids 418-518 of Human REV3L
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Store at -20°C. Stable for 12 months at -20°C

BACKGROUND

Introduction	Protein reversionless 3-like (REV3L) also known as DNA polymerase zeta catalytic subunit (POLZ) is an enzyme that in humans is encoded by the REV3L gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [5980](#)

Protein Refseq [NP_002903](#)

UniProt ID [O60673](#)
