



Anti-RGD1563533 (aa 1150-1199) polyclonal antibody (CABT-BL3144)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 1150-1199 (PQMPQIPSSV PHLPTSPLAT TSLESAKPQV KPGFLQFQDN DPCLATDCKY) of Rat RGD1563533 (XP_002726704).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [313713](#)