



Anti-RHOX11 (aa 69-110) polyclonal antibody (CABT-BL3152)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A 14 amino acid synthetic peptide found within the region of amino acids 69/110 of Mouse Rhox11 (VMPNTQDTGR EEPEETSKVA ETSEQSLFRI PRKAYRFTPG QLWELQAVFV)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	194738
Protein Refseq	NP_941000

UniProt ID

[Q810N8](#)
