



Anti-RIMS4 (internal region) polyclonal antibody (CABT-BL3154)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthesized peptide derived from an internal region of Human RIMS4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Cell junction; synapse.
Format	Liquid
Buffer	50% Glycerol, PBS, 0.88% Sodium chloride pH: 7.40
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

BACKGROUND

Introduction	Regulates synaptic membrane exocytosis.
--------------	---

GENE INFORMATION

Entrez Gene ID [140730](#)

Protein Refseq [NP_892015](#)

UniProt ID [Q9H426](#)
