



Anti-RNF220 (full length) polyclonal antibody (CABT-BL3166)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length human RNF220 protein (AAH98300.1)
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ICC/IF, WB
Format	Liquid
Size	50 µg
Buffer	1X PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	RING finger protein 220 is a protein that in humans is encoded by the RNF220 gene.
--------------	--

GENE INFORMATION

Entrez Gene ID [55182](#)

Protein Refseq [NP_060620](#)

UniProt ID [Q5VTB9](#)
