



Anti-RNF44 (aa 107-156) polyclonal antibody (CABT-BL3167)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the internal sequence amino acids 107-156 (LSYTVTTVT QGFPLPTGQH IPGCSAQQLP ACSVMFSGQH YPLCCLPPPL) of Human RNF44.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions.
---------------------	---

GENE INFORMATION

Entrez Gene ID	22838
Protein Refseq	NP_055716
UniProt ID	A0A024R7Q1