



# Anti-ROD1 (aa 81-130) polyclonal antibody (CABT-BL3174)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within N-terminal amino acids 81-130 (YYTPVTPHLR SQPVYIQYSN HRELKTDNLP NQARAQAALQ AVSAVQSGNL) of Mouse ROD1, isoform 2 (Q8BHD7).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Buffer</b>	98% PBS, 2% Sucrose
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	The protein encoded by this gene binds RNA and is a regulator of cell differentiation. The encoded protein;preferentially binds to poly(G) and poly(U) sequences in vitro. Several transcript variants encoding different;isoforms have been found for this gene.
---------------------	---

# GENE INFORMATION

Entrez Gene ID [230257](#)

---

Protein Refseq [NP\\_659153](#)

---

UniProt ID [G8JL74](#)

---