



Anti-ROPN1 (N-terminal) polyclonal antibody (CABT-BL3175)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	This antibody is specific for the N Terminus Region of the target protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC (PFA fixed)
Buffer	20 mM Potassium Phosphate, 150 mM Sodium Chloride, pH 7.0
Preservative	None
Storage	Store at -20°C. Avoid freeze-thaw cycles.

GENE INFORMATION

Entrez Gene ID	54763
Protein Refseq	NP_060048
UniProt ID	Q9HAT0