



Anti-RPAP2 (C-terminal) polyclonal antibody (CABT-BL3177)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from C-terminal region of Human RPAP2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus. Membrane.
Format	Liquid
Buffer	PBS, 0.88% Sodium chloride, 50% Glycerol pH: 7.40 Note: without Mg ²⁺ and Ca ²⁺
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID	79871
Protein Refseq	NP_079089

UniProt ID

[Q8IXW5](#)
