



Anti-RPD3 (aa 503-521) polyclonal antibody (CABT-BL3178)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: C-SGSGTGSGAGAKGAKENNI, corresponding to amino acids 503-521 of Fruit fly (<i>Drosophila melanogaster</i>) HDAC1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ICC/IF, ChIP, IP, WB
Cellular Localization	Nucleus.
Format	Liquid
Size	100 µl
Buffer	serum with 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [38565](#)

Protein Refseq [NP_647918](#)

UniProt ID [Q94517](#)
