



Anti-RPL39L (internal region) polyclonal antibody (CABT-BL3181)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to internal of human RPL39L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC, ELISA
Format	Liquid
Buffer	In PBS, pH 7.4 (150mM NaCl, 0.02% sodium azide, 50% glycerol)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	<p>This gene encodes a protein sharing high sequence similarity with ribosomal protein L39. Although the name of this gene has been referred to as ribosomal protein L39 in the public databases, its official name is ribosomal protein L39-like. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome.</p>
---------------------	--

GENE INFORMATION

Entrez Gene ID [116832](#)

Protein Refseq [NP_443201](#)

UniProt ID [Q96EH5](#)
