



Anti-RPN10 (full length) polyclonal antibody (CABT-BL3184)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length protein produced in E.coli.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Proteasome.
Format	Liquid
Size	100 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [856607](#)

Protein Refseq [NP_012070](#)

UniProt ID [P38886](#)
