



Anti-RPN12a (full length) polyclonal antibody (CABT-BL3185)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length protein produced in E.coli.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	842760
Protein Refseq	NP_176633
UniProt ID	Q9SGW3