



Anti-RPN13 (aa 11-110) polyclonal antibody (CABT-BL3186)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: RAGVCEYNEDSRLCTPIPVQGEIEIKPNEEEELGFWDFFEW RPTEKPVGRELDPISLILIPGETMWVPIKSSKSGRIFALV FSSNERYFFWLQEKNSGNLP, corresponding to amino acids 11-110 of <i>S. cerevisiae</i> Rpn13.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [851140](#)

Protein Refseq [NP_013525](#)

UniProt ID [O13563](#)
