



Anti-RPP25 polyclonal antibody (CABT-BL3188)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human RPP25.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC, IF
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Component of ribonuclease P, a protein complex that generates mature tRNA molecules by cleaving their 5-ends. Also a component of RNase MRP. This subunit binds to RNA.
--------------	---

GENE INFORMATION

Entrez Gene ID [54913](#)

Protein Refseq [NP_060263.2](#)

UniProt ID [Q9BUL9](#)
