



# Anti-RPP25 polyclonal antibody (CABT-BL3189)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugate	Unconjugated
Applications	WB, ELISA
Buffer	PBS with 0.09% (W/V) sodium azide
Preservative	0.09% Sodium Azide
Storage	RPP25 Antibody (C-term) can be refrigerated at 2-8°C for up to 6 months. For long term storage, place the at -2°C.

## BACKGROUND

Introduction	Component of ribonuclease P, a protein complex that generates mature tRNA molecules by cleaving their 5-ends. Also a component of RNase MRP. This subunit binds to RNA.
--------------	---

## GENE INFORMATION

Entrez Gene ID	<a href="#">54913</a>
----------------	-----------------------

Protein Refseq [NP\\_060263.2](#)

---

UniProt ID [Q9BUL9](#)

---