



Anti-RPS19BP1 (internal region) polyclonal antibody (CABT-BL3190)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from internal region of Human RPS19BP1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Nucleus; nucleolus.
Format	Liquid
Buffer	49% PBS, 0.88% Sodium chloride, 50% Glycerol pH: 7.40
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID	91582
Protein Refseq	NP_919307

UniProt ID

[A0A024R1T1](#)
