



Anti-RPS27L (aa 1-50) polyclonal antibody (CABT-BL3192)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MPLARDLLHP SLEEEKKKHK KKRLVQSPNS YFMDVKCPGC YKITTVFSHA) of Human RPS27L.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	The gene encodes a protein sharing 96% amino acid similarity with ribosomal protein S27,. This would suggest thst the encoded protein may be a component of the 40S ribosomal subunit.
---------------------	--

GENE INFORMATION

Entrez Gene ID [51065](#)

Protein Refseq [NP_057004](#)

UniProt ID [Q71UM5](#)
