



Anti-RPT1 (aa 1-100) polyclonal antibody (CABT-BL3194)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant (His-tagged) fusion protein: MPPKEDWEKY KAPLEDDDKK PDDDKIVPLT EGDIQVLKSY GAAPYAAKLK QTENDLKDIE ARIKEKAGVK ESDTGLAPSH LWDIMGDRQR LGEEHPLQVA, corresponding to amino acids 1-100 of Rpt1 (yeast).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. cerevisiae</i> , <i>S. pombe</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Size	50 µl
Buffer	Whole serum. Dilute to working strength with phosphate buffered saline pH 7.2-7.4 and 1% normal goat serum (if a goat anti-rabbit IgG linker antibody is to be used).
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [853712](#)

Protein Refseq [NP_012777](#)

UniProt ID [P33299](#)
